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1. Oral, E, Muratoglu, O, "The effects of high dose irradiation on the cross-linking of vitamin E-blended ultrahigh molecular weight polyethylene," *Biomaterials* 29 (2008) 3557-3560
2. Oral, E, Muratoglu, O, "Vitamin E diffused, highly crosslinked UHMWPE: a review," *Int Orthop*. 2011 Feb; 35(2): 215-223
3. Exactechの社内資料
4. Muratoglu OK, Bragdon CR, O'Connor DO, Perinchief RS, Estok DM, Jasty M, Harris WH (2001) Larger diameter femoral heads used in conjunction with a highly cross-linked ultra-high molecular weight polyethylene: A new concept. *J Arthroplasty* 16(8 Suppl):24 -30
5. Oral E, Christensen S, Malhi A, Wannomae K, Muratoglu O (2006) Wear resistance and mechanical properties of highly crosslinked UHMWPE doped with vitamin E. *J Arthroplasty* 21 (4):580 -591
6. Muratoglu OK, Bragdon CR, O'Connor DO, Jasty M, Harris WH, Gul R, McGarry F (1999) Unified wear model for highly crosslinked ultra-high molecular weight polyethylenes (UHMWPE). *Biomaterials* 20(16):1463 -1470
7. Ries M, Pruitt L (2005) Effect of crosslinking on the microstructure and mechanical properties of ultra-high molecular weight polyethylene. *Clin Orthop Relat Res* 440:149 -156
8. Oral E, Wannomae KK, Hawkins NE, Harris WH, Muratoglu OK (2004) a-Tocopherol doped irradiated UHMWPE for high fatigue resistance and low wear. *Biomaterials* 25(24):5515 -5522
9. Bhattacharyya, S, Spiegelberg, S.H, Harris, W.H, Muratoglu, O.K. "Mechanical Elimination of Free Radicals in UHMWPE After Radiation." 50th Annual Meeting of the Orthopaedic Research Society. Paper No: 0184
10. Ebru Oral, Bassem W. Ghali, Orhun K. Muratoglu. "The elimination of free radicals in irradiated UHMWPEs with and without vitamin E stabilization by annealing under pressure." *J Biomed Mater Res B Appl Biomater*. 2011 Apr;97(1):167-74. doi: 10.1002/jbm.b.31799. Epub 2011 Feb 24

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本書で言及される製品は、国によって異なる商標で販売されている可能性があります。すべての著作権、申請中および登録済みの商標は、Exactechに帰属します。本資料は、Exactech販売担当者および医師の使用および両者に利することのみを意図しています。Exactechによる明示的な書面許可なしに本書を再配布、複写、開示できません。©2019 Exactech. 12-0000586 0819